

## LifeLike Long Femoral Leg Simulator Set-Up Instructions

### Parts Included

- Line Tubing x 3
- Luer Lock Vessel Connectors A (larger) x 2
- Luer Lock Vessel Connectors B (smaller) x 4
- Reservoir with Push-to-Connect connectors
- Pump with Push-to-Connect connectors
- Zip ties
- Extra tubing and connectors

### Assembly

1. Fill the reservoir with water and screw on the lid.
2. Attach the **blue** end of the short piece of tubing to one of the **blue** connectors on the reservoir by pushing the tubing into the connector. Push the tubing in so the grey piece of plastic in the tubing is all the way into the connector.
3. Attach the **green** end of the same tube to the **green** inlet connector on the pump. If looking at the face of the pump, the inlet connector is on the left.
4. Attach the **green** end of one of the long pieces of tubing to the **green** outlet connector on the pump. If looking at the face of the pump, the inlet connector is on the right.
5. Attach the **red** end of the same tube to the **red** connector in the upper thigh area of the simulator.
6. Attach the **black** end of the remaining piece of tubing to the **black** connector in the ankle area of the simulator.
7. Attach the **blue** end of the same tube to the open **blue** connector on the reservoir.
8. Remove the LifeLike Medium Femoral Artery model from its packaging. Insert a **Luer Lock Vessel Connector A** into one end of the vessel and secure the connector with a zip tie. Place the vessel into the Groin cavity and roughly measure the length needed to reach the connector at the other end. Trim the vessel to size, then insert a **Luer Lock Vessel Connector A** into the other end of the vessel and secure the connector with a zip tie. Repeat for the Small Femoral Artery model using **Luer Lock Vessel Connector B** for each. (3 models needed total per training session: [1] Medium Femoral and [2] Small Femorals).
1. To insert the Femoral Arteries into the openings in the leg simulator, line up the Luer Lock Vessel Connector with the Inner Cavity Luer Lock Connector and give a quarter turn clockwise to connect. For the other end, twist the vessel counter-clockwise a quarter turn, insert into the connector on the opposite side of the leg opening, then turn clockwise to connect. This should allow the vessel to be inserted without any twisting or kinking of the vessel. If the vessel is twisted, remove and try again. Repeat for the remaining cavities.
9. Plug in the pump, turn on, and let the water run through the system while checking for any leaks.


**\*NOTE\*:** Extra connectors have been provided to set up the vessels prior to the previous training session being completed for a quick switch out. Simply insert the connectors and secure with zip ties. When the training session is completed, remove the used vessels from the simulator, cut away the zip ties and remove the connectors. The used vessels can be discarded and prepared vessels can be inserted.

### Disassembly

1. Turn off the pump. Loosen the lid on the reservoir slightly. Remove the tube attached to the pump from the reservoir by holding down the grey ring on the outside of the connector while pulling the tubing outwards. Hold tube end in the air so water does not spill out.
2. Turn on the pump and run until water has emptied from the entire system back into the reservoir. Once empty, turn off the pump and unplug.
3. Disconnect the used vessels from the simulator. Cut away the zip ties and remove the connectors.
4. Disconnect all tubing and empty the reservoir.
5. Dry off the simulator before storing.

### Diagram

## LIFELIKE LONG FEMORAL LEG SIMULATOR SET-UP



#### LEGEND

- A - OUT FROM RESERVOIR
- B - INTO PUMP
- C - OUT FROM PUMP
- D - INTO LEG SIMULATOR
- E - OUT OF LEG SIMULATOR
- F - INTO RESERVOIR